

9)))))

Innovation with engineering legacy



Artemis Proximal Femoral Nail System

The Artemis Proximal Femoral Nail System is indicated for fixation of stable and unstable intertrochanteric fractures and subtrochanteric, including but not limited to nonunion, malunion and tumor resections.



Oart



Design Features

The Artemis Proximal Femoral Nail System encompasses four implant components. All components are available in sterile packaging only.

Proximal Femoral Nail Kit

The proximal femoral nail kits, consist of the **Nail** and a preassembled **Set Screw**. Its universal design is intended for left and right application.

The nail is constructed of a titanium alloy core encompassed by injection molded carbon fiber reinforced (CFR) polyether ether ketone (PEEK). The set screw is made of titanium alloy.

Lag Screw

The lag screw has a diameter of 11mm and is available in lengths ranging from 70 to 130mm in 5mm increments. The most common lengths are offered in 2.5mm increments. The lag screw is made of titanium alloy.

Anti-Rotational Locking Pin


The anti-rotational locking pin is designed to provide additional rotational stability and fixation performance.

Locking Screw

The locking screw of the Artemis Proximal Femoral Nail System is blunt tipped, self tapping, dual lead thread and have a diameter of 5mm. It is fully threaded and comes in a wide range from 25mm to 110mm in 2.5mm increments up to 50mm and in 5mm increments up to 110mm. A high torque transfer is achieved through the T25 torx screw head.


Instrumentation

The Artemis system has a state-of-the-art instrument platform. The instruments are designed for a minimally invasive surgical technique.




Anti-Rotational Locking Pin (70, 95, 110mm) Ø5mm

35mm



Locking screw



300-440mm

increments)

(in 20mm



Length dependent RoC

Length (mm)	Radius (m)
300	1.00
320	1.05
340	1.10
360	1.15
380	1.20
400	1.25
420	1.30
440	1.35







Added Visibility

The composite materials provide clinical benefits with improved visualization in both long and short nail configurations. This added visibility allows for easier assessment of bone healing and Lag Screw position.

Hybrid Materials and Manufacturing

Patented* manufacturing process utilizes a combination of titanium and composite materials. The nails structural benefits are provided by the strength of titanium and the flexibility of CFR PEEK which allows for increased micro movements and callus formation. The process of minimized titanium milling with the benefit of injection molded CFR PEEK results in reduced waste manufacturing. Titanium

CFR PEEK



Lagshield[™] Function

Lagshield is a protective layer of CFR PEEK on the lateral side of the nail.

This prevents the lag screw reamer from notching the nails load bearing titanium core during the reaming procedure.




The proximal femoral nail kits, consist of the **Nail** and a preassembled **Set Screw**. Its universal design is intended for left and right application.



Smooth nail insertion, preventing anterior impingement

Lagshield[™] Function

Lagshield is a protective layer of CFR PEEK on the lateral side of the nail. This prevents the lag screw reamer from notching the nails load bearing titanium core during the reaming procedure.



Lag shield aids in guidance of the reamer and protects the load bearing titanium, preventing the lag screw reamer from notching the nails titanium core during the reaming procedure.



Stability and Rotational Control

Enhanced Stability and Rotational Control: Consists of a Lag Screw with optional combination of integrated anti-rotational Locking Pin which supports compression and helps to remove the Z-effect.



A simple, one piece cannulated set screw locks in to the lag screw, preventing rotation of the lag screw.



The lag screw can move backwards, with the head fragment but cannot rotate.



Locking pin locks into the lateral cortex and the nail, providing additional stability allowing the head fragment to collapse with the lag screw.





Added Visibility The composite materials provide clinical benefits with improved visualization in both long and short nail configurations. This added visibility allows for easier assessment of bone healing.

The radiolucent proximal end of the nail provides better radiographic visibility of the femoral neck in an oblique medio-lateral image.









This material is intended for health care professionals. Distribution to any other recipient is prohibited. For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert.

This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, A written operative technique is available at Innov8ortho website or through your local Innov8ortho representative. As with all surgical procedures, the technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure. Caution: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only.

GLW, Inc, GLW Medical Innovation, and Artemis are trademarks of GLW, Inc. Innov8ortho, LLC is the exclusive distributor of the Proximal Femoral Nail System. Patents: www.glwmed.com/patents.

For additional information on Innov8ortho® and its products please visit: www.innov8ortho.com



Scan for more product information



Distributed by:

Innov8ortho, LLC 300 Sylvan Ave, 2nd Floor Englewood Cliffs, NJ 07632

custsvc@innov8ortho.com



GLW, Inc. 300 Sylvan Ave, 2nd Floor Englewood Cliffs, NJ 07632