



Preparation for Nail Insertion





● Non Sterile Packaged
● Sterile Packaged

Insertion of the Nail











Preparation for Lag Screw Insertion

23	24	25	26	27	28	29	30	31	32	33	34	35
												
Top Inset Tray T5-0110-000 T-Handle, Zimmer/Hall Coupling	Middle Inset Tray T5-0240-000 K-wire Sleeve	Middle Inset Tray T5-0230-000 Lag Screw Guide Sleeve	Procedure Kit Short Nail T4-0200-010S K-wire, Threaded Trocar Tip, 3.5 x 450mm	Middle Inset Tray 5-0280-000 Lag Screw Measuring Gauge	Middle Inset Tray T5-0270-000S Anti-Rotation Wire	Top Inset Tray T5-0100-000 Handle, AO small	Base Inset T5-0390-000 Screwdriver Bit T25, AO Small	Middle Inset Tray T5-0250-000 Lag Screw Drill, Ø11mm, Zimmer/Hall	Middle Inset Tray T5-0260-000 Lag Screw Drill Stop	Middle Inset Tray T5-0290-000 Lag Screw Driver	Middle Inset Tray T5-0300-000 Fixation Rod for Lag Screw Driver	Base Inset T5-0530-000 Lag Screw Tap

● Non Sterile Packaged
● Sterile Packaged

Insertion of the Lag Screw

36	37	38	39	40	41	42	43
							
Middle Inset Tray T5-0290-000 Lag Screw Driver	Middle Inset Tray T5-0300-000 Fixation Rod for Lag Screw Driver	Sterile Packaged T6-0200-0XX Lag Screw	Base Inset T5-0390-000 Screwdriver Bit T25, AO Small	Top Inset Tray T5-0100-000 Handle, AO small	Middle Inset Tray T5-0450-000 Compression Device	Base Inset T5-0310-000 Flexible Screwdriver, 4.0mm, AO Small	Top Inset Tray T5-0100-000 Handle, AO small

Preparation for Locking Screw Insertion

Insertion of the Locking Screw

44




Base Inset
T5-0340-000
Locking
Screw Trocar

45



Base Inset
T5-0330-000
Locking Screw
Drill Sleeve

46



Base Inset
T5-0320-000
Locking Screw
Guide Sleeve

47



**Procedure Kit
Short Nail**
T4-0200-010S
Locking Screw Drill,
Long

48



Top Inset Tray
T5-0100-000
Handle, AO small

49



Base Inset
T5-0390-000
Screwdriver Bit T25,
AO Small

50



**Sterile
Packaged**
T6-0300-0XX
Locking Screw

51



Base Inset
T5-0220-000
Ball Tip Screwdriver,
8.0mm, Zimmer/Hall

52



Top Inset Tray
T5-0110-000
T-Handle, Zimmer/
Hall Coupling

● Non Sterile Packaged
● Sterile Packaged

Extraction

53



Top Inset Tray
T5-0100-000
Handle, AO small

54



Base Inset
T5-0390-000
Screwdriver Bit T25,
AO Small

55



Top Inset Tray
T5-0290-000
Lag Screw Driver

56



Middle Inset Tray
T5-0300-000
Fixation Rod for
Lag Screw Driver

57



Top Inset Tray
T5-0100-000
Handle, AO small

58



Base Inset
T5-0310-000
Flexible Screwdriver,
4.0mm, AO Small

59



Base Inset
T5-0430-000
Extraction Rod,
Zimmer/Hall

60



Base Inset
T5-0440-000
Slotted Hammer

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a surgeon. Rx only.

The information presented is intended to demonstrate a Innov8ortho product. A surgeon must always refer to the package Inset, product label and/or instructions for use, including the instructions for Cleaning and Sterilization (if applicable), before using any Innov8ortho product.

GLW, Inc, GLW Medical Innovation, and Artemis are trademarks of GLW, Inc. Innov8ortho, LLC. is the exclusive distributor of the Artemis Proximal Femoral Nail System.

Distributed by:

Innov8ortho, LLC.
300 Sylvan Ave, 2nd Floor
Englewood Cliffs, NJ 07632